



FILTERS, SEGMENTERS, GRIDS & MATTRESSES

FILTERS

<i>Code</i>	<i>Description</i>	<i>Units</i>
<i>C06087C</i>	<i>FILTER DECUBITUS (VARIOUS AREAS)</i>	<i>EACH</i>
<i>C06087D</i>	<i>FILTER HIP TO ANKLE (SWIMMERS NECK)</i>	<i>EACH</i>
<i>C06087E</i>	<i>FILTER WHOLE SPINE AP (LATERAL THORACIC LUMBER SPIN)</i>	<i>EACH</i>
<i>C06087</i>	<i>FILTER THORACIC SPINE</i>	<i>EACH</i>
<i>C06087F</i>	<i>FILTER SHOULDER</i>	<i>EACH</i>
<i>C06087B</i>	<i>FILTER FOOT</i>	<i>EACH</i>
<i>C06087G</i>	<i>FILTER HIP TO ANKLE</i>	<i>EACH</i>

SEGMENTERS

<i>Code</i>	<i>Description</i>	<i>Units</i>
<i>C06071</i>	<i>SEGMENTER 18 x 24cm Long Axis, Short Axis, or L-Shape</i>	<i>EACH</i>
<i>C06074</i>	<i>SEGMENTER 24 x 30cm Long Axis, Short Axis, or L-Shape</i>	<i>EACH</i>
<i>C06073</i>	<i>SEGMENTER 35 x 43cm Long Axis, Short Axis, or L-Shape</i>	<i>EACH</i>
<i>C06072</i>	<i>SEGMENTERS QUARTERS, CUT OUT CIRCLES, FILM FRAMING</i>	
<i>C06073A</i>	<i>POLYCARBONATE BASE/VELCRO</i>	
<i>C06073B</i>	<i>POLYCARBONATE L-SHAPED PATIENT ARM BOARDS (ANGIO CT)</i>	

GRIDS

<i>Code</i>	<i>Description</i>	<i>Units</i>
<i>EX9011</i>	<i>GRIDS 85 / 103 LINES 8:1 RATIO (18x24, 24x30, 30x40, 35x35, 35x43)</i>	<i>EACH</i>
<i>EX9012</i>	<i>GRIDS 85 / 103 LINES 10:1 RATIO (ALL SIZES AVAILABLE)</i>	<i>EACH</i>
<i>EX9014</i>	<i>GRIDS 85 / 103 LINES 12:1 RATIO (ALL SIZES AVAILABLE)</i>	<i>EACH</i>
<i>EX9013A</i>	<i>GRIDS PERSPEX COVERS</i> <i>Short, Medium & Long Focus available</i> <i>OTHER SIZES ARE AVAILABLE UPON REQUEST</i>	

MATTRESSES

<i>Code</i>	<i>Description</i>	<i>Length x Width x Depth</i>	<i>Units</i>
<i>C06021</i>	<i>MATTRESS</i>	<i>1840 x 450 x 25mm</i>	<i>EACH</i>
<i>C06022</i>	<i>MATTRESS</i>	<i>1840 x 450 x 50mm</i>	<i>EACH</i>
<i>C06023</i>	<i>MATTRESS</i>	<i>2140 x 450 x 25mm</i>	<i>EACH</i>
<i>C06023A</i>	<i>MATTRESS</i>	<i>2140 x 600 x 25mm</i>	<i>EACH</i>
<i>C06024</i>	<i>MATTRESS</i>	<i>1840 x 600 x 25mm</i>	<i>EACH</i>
<i>C06024A</i>	<i>MATTRESS</i>	<i>2000 x 650 x 25mm</i>	<i>EACH</i>
<i>C06024B</i>	<i>MATTRESS</i>	<i>1800 x 600 x 25mm</i>	<i>EACH</i>
<i>C06024C</i>	<i>MATTRESS</i>	<i>2400 x 750 x 25mm</i>	<i>EACH</i>
<i>C06024D</i>	<i>MATTRESS</i>	<i>1840 x 600 x 50mm</i>	<i>EACH</i>

**ALL MATTRESSES ARE; VINYL COVERED, ANTI-STATIC, ANTI-BACTERIAL, ANTI FUNGICIDAL
ALTERNATE SIZES ARE AVAILABLE UPON REQUEST**